EQ-1 Handheld Reader  
Lecteur portable EQ-1  
Lector portátil EQ-1  
Aparelho portátil EQ-1  
Lettore portatile EQ-1  
EQ-1 Handlesegerät

# USER'S MANUAL

MANUEL D'UTILISATION  
MANUAL DEL USUARIO  
MANUAL DO USUÁRIO

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# Contents

---

|  |    |
|--|----|
| <b>1. Stablelab EQ-1 Reader Box Contents</b> .....   | 4  |
| <b>2. System Components</b> .....                    | 5  |
| <b>3. Getting Started</b> .....                      | 6  |
| <b>4. Sample Preparation</b> .....                   | 7  |
| 4.1 Whole blood samples                              |    |
| 4.2 Serum and Plasma samples                         |    |
| <b>5. Standard Mode</b> .....                        | 9  |
| <b>6. Quick Mode</b> .....                           | 10 |
| <b>7. Result Overview</b> .....                      | 11 |
| <b>8. Maintenance</b> .....                          | 12 |
| 8.1 How to change Reader Battery                     |    |
| 8.2 "Clean Lens" error: How to clean the Reader Lens |    |
| <b>9. Error Message/Troubleshooting</b> .....        | 13 |
| 9.1 Battery Low                                      |    |
| 9.2 Clean Lens                                       |    |
| 9.3 Check Insert                                     |    |
| 9.4 Too Late   |    |
| 9.5 No C-Line  |    |
| 9.6 Calibrate Pass                                   |    |
| 9.7 Calibrate Fail                                   |    |
| <b>10. Precautions</b> .....                         | 14 |
| 10.1 Intended Use                                    |    |
| 10.2 Reader Care                                     |    |
| 10.3 Operating/Storage Conditions                    |    |
| 10.4 Disclaimer                                      |    |
| 10.5 Limitations of Use                              |    |
| 10.6 Symbol Descriptions                             |    |
| <b>11. Support</b> .....                             | 16 |

**FOR VETERINARY USE ONLY**

## 1. Stablelab EQ-1 Reader Box Contents

## 2. System Components


Front View

Bottom View



Side View

- 1 Reader Display
- 2 Reader Mode Switch
- 3 Cartridge Port Door
- 4 Cartridge Port
- 5 Battery Port
- 6 Battery Door

### 3. Getting Started


1. Insert the Duracell® Coppertop 9V battery provided, ensure the + and - symbols on the battery and EQ-1 Reader align (a).
2. Use the enclosed screwdriver to attach and close the battery door (b).  
**TIP:** The ridge on the battery door should lay directly over the + and - symbols and align with the ridge in the middle of the side of the reader (c).
3. Move the reader mode switch to either Standard or Quick Mode (refer to sections 5 & 6).

## 4. Sample Preparation

The volume of sample used is critical for whole blood, serum and plasma samples. For best results, use fresh whole blood immediately or refrigerated whole blood, serum and plasma samples within 12 hours of collection.


### 4.1 Whole blood samples

- Remove the Mix Solution bottle cap. Firmly insert the BloodCaptor collection arm with blood-filled bulb completely into the Mix Solution so that a tight seal is formed (b).
- Gently agitate the contents and run the sample immediately according to the instructions on the Stablelab Serum Amyloid A Test package insert (c).

- Whole blood samples are the easiest, quickest and most convenient to use at the horse's side.
- Use standard venipuncture with appropriately sized needle and syringe to obtain sample.
- Utilize the BloodCaptor™ to collect the sample directly from the syringe hub via capillary action, immediately following collection (preferred) (a).
- Alternatively, whole blood can be used from any anticoagulated blood collection tube.
- A BloodCaptor is included with each individual test pack and accurately measures 8.5  $\mu$ l of whole blood.
- Ensure the BloodCaptor bulb is full but avoid excess blood on the outside of the bulb. Should this occur, carefully wipe off, taking care not to disturb the collected sample within the BloodCaptor.

## 4.2 Serum and Plasma samples

- Any standard blood collection tube can be used including lithium heparin (green top), EDTA (purple top) and sodium citrate (blue tops) tubes, as well as serum tubes (red/tiger top).
- After appropriate centrifugation, use a micropipette device that can measure 5  $\mu$ l, in place of the BloodCaptor.
- Measure 5  $\mu$ l of serum or plasma from the blood tube and deposit into the Mix Solution bottle (d).
- Firmly insert the empty BloodCaptor collection arm completely into the Mix Solution so that a tight seal is formed (e).
- Gently agitate the contents and run the sample immediately according to the instructions on the Stablelab Serum Amyloid A Test package insert (f).

## 5. Standard Mode

Standard Mode is the simplest and quickest way to run one Stablelab SAA cartridge in 10 minutes. Your EQ-1 Reader will time the test and show the result on the screen after 10 minutes.

a

b

1. Prepare the sample by following the instructions on the Stablelab Serum Amyloid A Test package insert.
  - Test cartridges should have sample applied within 10 minutes of opening the foil packaging.
2. Ensure the reader mode switch is pushed to the left (a).
3. Insert a new Stablelab SAA test cartridge into the cartridge port of the reader with the result window of the cartridge facing upwards and the sample well protruding from the cartridge port.
4. The EQ-1 Reader will automatically turn on and display “Apply Sample”.
5. To apply sample to the cartridge, gently squeeze the bottle and discard the first few drops. Then apply exactly 4 drops from the Mix Solution bottle to the cartridge well.
6. The reader will detect the sample within 60 seconds of sample application and will initiate an automatic 9-minute countdown (b).
  - Until the countdown is initiated, “Apply Sample” will remain visible on the screen - do not apply additional drops.
7. After a total of 10 minutes, the result is displayed for 20 seconds before reader powers down to conserve energy.
8. The result is stored until the cartridge is removed. The result will appear for 5 seconds after cartridge removal before being deleted from the device.
9. Reader will power down on its own to conserve energy.

## 6. Quick Mode

Quick Mode is used to run multiple different Stablelab SAA cartridges in sequence. Each test is timed by you independent of the reader, requiring precise 10-minute timing/cartridge. The cartridges are individually inserted in series into the reader at the end of 10 minutes per cartridge to provide a quantitative SAA value.

a

b

1. Prepare the sample by following the instructions on the Stablelab Serum Amyloid A Test package insert.
  - Test cartridges should have sample applied within 10 minutes of opening the foil packaging.
2. Ensure the reader mode switch is to the right. The letter Q will be displayed in the upper right corner of the screen when the cartridge is inserted (a).
3. Have one Stablelab SAA test cartridge per sample available.


4. When running multiple samples, label the Mix Solution bottle and respective test cartridge to avoid confusion.
5. To apply sample to the cartridge, first gently squeeze the bottle and discard the first few drops. Then apply exactly 4 drops from the bottle to the cartridge well. Set a timer for 10 minutes. Repeat procedure for each additional sample (b).
6. As the sample moves up the test cartridge a visible control line at 'C' will appear, indicating that the test is working.
7. Keep cartridges on a flat horizontal surface while developing. Exactly 10 minutes after sample application, insert the cartridge into the cartridge port with the result window of the cartridge facing upwards and the sample well protruding from the cartridge port.
8. SAA result will be displayed immediately and stays on screen for 20 seconds before being deleted permanently from the device.
9. After the SAA result is displayed and recorded, remove the test cartridge and insert the next as each additional test cartridge reaches the 10-minute mark.
10. **It is critical to read each test cartridge at exactly 10 minutes** because the cartridge Control and Test lines continue to develop over time. Inserting a cartridge beyond the 10-minute mark will falsely elevate the SAA value measured. Conversely, inserting a cartridge before 10 minutes have passed, will read a lower SAA value than is present.
11. Reader will power down on its own to conserve energy.

## 7. Result Overview

The reader will display the result in  $\mu\text{g/ml}$  on the screen after 10 minutes (in Standard Mode) or after the test cartridge has been inserted after exactly 10 minutes have passed (in Quick Mode).

The detection range is between 0 - 3000  $\mu\text{g/ml}$ , where zero means SAA is not detected.


- Test cartridges **must** be read by the EQ-1 reader at exactly 10 minutes because the cartridge Control and Test lines continue to develop over time.
- This is especially important in Quick Mode as the test is timed by you, independent of the reader.
- Inserting a cartridge beyond the 10-minute mark will display a falsely elevated SAA value.
- Conversely, inserting a cartridge before 10 minutes have passed, will read a lower SAA value than is present.
- **Important:** Re-inserting the test cartridge beyond the 10-minute mark, to recheck the result is not recommended and will give an erroneous, falsely elevated result.

## 8. Maintenance

If your EQ-1 Reader is not performing optimally, reset the device. To reset the device, simply remove the battery and then reinsert or replace (see section 8.1 for more information)

### 8.1 How to change the Reader Battery

1. The battery port is located on the left side of the Reader.
2. Remove the two screws from the battery door using the screwdriver provided (a).
3. Pay close attention to replace the old battery with a new 9V battery, taking care to align the + and - symbols with the battery port (b).
  - For optimal performance, a Duracell Coppertop 9V disposable battery is recommended.
4. Replace the battery door using the screwdriver provided, ensuring that it has been inserted in the correct orientation.
  - **Important:** The ridge on the battery door should lay directly over the + and - symbols and align with the ridge in the middle of the side of the reader (c).
  - Do not overtighten the screws, as they may break or strip.

### 8.2 “Clean Lens” error: How to clean the Reader Lens

1. Open the cartridge port door.
2. Tilt the reader up to locate the lens through the cartridge port door on the inside upper surface of the Reader.
3. The lens is a silver-black strip coursing front to back in the middle of the white rectangle.
4. Using a clean cotton swab dipped in 70% isopropyl alcohol, carefully clean the lens by running the cotton swab across the lens in a forward and backwards motion. Do not apply any other liquid or cloth types directly to the lens.
5. After the lens has been cleaned, examine to ensure no particulates remain.
6. **Important: As the final step, the reader MUST be Reset.** To do this, remove the battery from the device completely, re-insert or replace it (see section 8.1).
7. The reader is now ready for use.

## 9. Error Message/Troubleshooting

| Display Message           | What it means   | What to do   |
|---------------------------|---|--|
| <b>9.1 Battery Low</b>    | Indicates that your EQ-1 Reader is running low on battery and should be changed immediately.  | Follow instructions in section 8.1 for more information.   |
| <b>9.2 Clean Lens</b>     | The reader has detected dirt on the lens that may be affecting performance.   | Follow instructions in section 8.2 for more information.   |
| <b>9.3 Check Insert</b>   | Appears when a Stablelab SAA test cartridge has not been inserted correctly into the Reader.  | Ensure that the Stablelab SAA cartridge is pushed fully into the cartridge port with the result window facing upwards and the sample well protruding from the cartridge port.  |
| <b>9.4 Too Late</b>       | Appears in Standard Mode when your EQ-1 Reader detects that a used Stablelab SAA test cartridge has been inserted.  | To read a used Stablelab SAA test cartridge, switch the EQ-1 Reader to Quick Mode (see section 6 for more information). If a new Stablelab SAA cartridge has been inserted in standard mode and the message is displayed, please contact Zoetis for further technical support. |
| <b>9.5 No C-Line</b>      | A "No C-Line" message indicates the test cartridge has failed to run. This occurs by either too little or too much sample being applied to the test cartridge well or sample being applied in operating temperatures outside recommendations. | Prepare the sample by following the instructions on the Stablelab Serum Amyloid A Test package insert and sections 4-6. If, on examination of a test cartridge, a C- Line is present, please contact Zoetis for further technical support.                                     |
| <b>9.6 Calibrate Pass</b> | A new, blank SAA test cartridge has been inserted in Quick Mode initiating a calibration check and the check has passed.  | Continue to use your EQ-1 Reader as normal and prepare the sample by following instructions on the Stablelab Serum Amyloid A Test package insert.  |
| <b>9.7 Calibrate Fail</b> | A new, blank SAA test cartridge has been inserted in Quick Mode initiating a calibration check and the check has not passed.  | Contact Zoetis for further technical support.  |

## 10. Precautions

### 10.1 Intended Use

This device provides an *in-vitro* quantitative determination of equine Serum Amyloid A (SAA). This device is intended for use with equine whole blood, serum or plasma using Stablelab SAA test cartridges only. This device is an analytical, point of care device designed for veterinary use only.

### 10.2 Reader Care

- Do not operate close to fire or after the device has contacted water.
- Do not expose or immerse the device in liquid.
- The EQ-1 Reader is only for use with Stablelab SAA cartridges. Inserting incompatible products may damage the device.
- Do not allow foreign objects to enter the device.
- Do not use if the device is damaged.
- Do not attempt to disassemble or repair the reader. Failure to observe this precaution could result in fire or electrical shock. If the device needs repair, contact your local Zoetis representative and/or distributor.

### 10.3 Operating/Storage Conditions

- This device is calibrated to operate between 15°C - 35°C. Using the device outside of this temperature range may affect accuracy. Do not operate the device in direct sunlight.
- Performance of the EQ-1 Reader may be impaired if the device is stored below -10°C or above 50°C. Store the device in a cool, dry place. Do not store in locations exposed to extreme temperatures or moisture.

### 10.4 Disclaimer

No part of this publication may be reproduced or transmitted in any form or by any means, electronic or mechanical, including photocopy, recording, or any information storage or retrieval system, without permission in writing from Zoetis, Inc. Zoetis, Inc. reserves the right to modify its products and services at any time. Information in this User's Manual is subject to change without notice. Although prepared to ensure accuracy, Zoetis assumes no liability for errors, or for any damages resulting from application or use of this information. No rights are granted for diagnostic uses. No rights are granted for use of the product for commercial application of any kind, including but not limited to manufacturing, quality control, or commercial services, such as contract services or fee for services. It is the responsibility of the purchase/end user to acquire any additional intellectual property rights that may be required.

### 10.5 Limitations of Use

Only Stablelab SAA test cartridges are to be used with the device. The device has a detection range of 0-3000 µg/ml of equine SAA. Factors that might affect performance of the test include the quality and volume of the sample, accuracy of the reagents, sample handling, time, and temperature. Great care has been taken to ensure the quality and the reliability of this product. However, it is possible that in certain cases unusual results may be obtained due to a high level of interfering factors.

## 10.6 Symbol Descriptions

| Symbol  | Description                  | Symbol  | Description   |
|---|------------------------------|---|---|
|  | Keep away from sunlight.     |  | Manufacturer  |
|  | Temperature limitations.     |  | Expiration Date                                     |
|  | Consult instructions for use |  | Batch code (Lot)                                    |
|  | Serial number                |  | Authorized Representative in the European Community |
|  | WEEE Directive 2002/96/EC    |  | Denotes conformity to specific European directives  |

## 11. Support

Sales Order & Technical support

For United States, please contact your Zoetis Representative, Distribution Partner or Zoetis directly at T 1-888-963-8471

If you are based in Australia, Brazil, Canada, France, Netherlands or the United Kingdom to order please contact Zoetis directly using the customer service details below:

Australia:

[customerservice.au@zoetis.com](mailto:customerservice.au@zoetis.com)

T +1800 022 442

Brazil

[adm-sac@zoetis.com](mailto:adm-sac@zoetis.com)

T 0800 011 1919

Canada

[StablelabCanada@zoetis.com](mailto:StablelabCanada@zoetis.com)

T 1-800-663-8888

France

[commandes.zoetis@zoetis.com](mailto:commandes.zoetis@zoetis.com)

T +0810 734 937

Netherlands

[NLorderdesk@zoetis.com](mailto:NLorderdesk@zoetis.com)

T + 0845 300 8034

Germany

[bestellannahme@zoetis.com](mailto:bestellannahme@zoetis.com)

T + (0) 30 33 00 63 - 331

United Kingdom

[Zoetis-DirectAccountUK@zoetis.com](mailto:Zoetis-DirectAccountUK@zoetis.com)

T + 0845 300 8034

For all other countries, please contact your local market Zoetis representative.

Learn more about Stablelab by visiting [www.stablelab.com](http://www.stablelab.com)

Zoetis Inc., 333 Portage Street, Kalamazoo, MI, 49007, USA

Zoetis Belgium, Rue Laid Burniat 1, 1348 Louvain-la-Neuve, Belgium

**[www.stablelab.com](http://www.stablelab.com)**

All trademarks are the property of Zoetis Services LLC or a related company or a licensor unless otherwise noted. Stablelab is a registered trademark of Epona Biotech Limited, used under license. Duracell is a registered trademark of Duracell, Inc. © 2020 Zoetis Services LLC. All rights reserved.

The Zoetis logo is rendered in a bold, orange, lowercase, sans-serif font. The letter 'z' is stylized with a thick, curved underline that extends under the 'o'.

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Zoetis Inc., 333 Portage Street, Kalamazoo, MI, 49007, Stati Uniti