

**It's not magic. It's magnetic.**



# AlteCap™ MagCap™ System

TRANSFORMING THE PROCESS OF DECAPPING AND RECAPPING SAMPLES



## Making History

For the first time, equipment engineers and plastic experts have collaborated to design decappers and sample storage tubes in tandem.

## Made for each other

Integrating magnets into cap drivers eliminates the necessity for friction to secure a cap, unlike conventional decapping systems.

## Increasing Productivity

Dependable decapping and recapping cycles without issues related to cap wear and tear, improperly sealed tubes, or caps dislodging.

