

NEW AlteCap™ Switch

Semi-automated screw cap handling for any manufacturer's SBS tube types

Quickly cap, decap or recap tubes by column or row. Featuring interchangeable cap driver cassettes for easy processing of different screw cap tube types and tubes in 24, 48 or 96 SBS rack formats.

Part No.	Product name
----------	--------------

404000	AlteCap Switch
--------	----------------



Unparalleled design

Innovative two-handed control, for maximum user comfort and effortless screw cap processing

Intuitive to use

Light touch buttons initiate decapping or capping. Low force dual button cap ejection

Uniquely flexible

Interchangeable cap driver cassettes compatible with screw cap tubes in 24, 48 or 96 SBS rack formats from any manufacturer

Functionality

Process caps in columns or rows. Recap tubes or discard used caps. Apply fresh caps from a cap carrier

Precise torque control

Individual internal cap driver controls ensure each screw cap is perfectly capped every time

Multi use stand

Includes portrait or landscape rack placement. Convenient for storage, USB charging and holding caps for re-capping

AlteCap™ Switch

Semi-automated handheld decapper

The only decapper capable of switching between 6, 8 and 12 formats

Traditional handheld decappers are uncomfortable to hold, restrict the view of tubes during processing and only work with one tube and rack type. The revolutionary designed AlteCap™ Switch overcomes all of these limitations.

Balanced
two-handed
controls

Comfortable
ergonomic
curved design

Clear view
of caps while
processing



Switch cassettes for
different tube or
rack types from any
manufacturer

Perfect
capping of
every tube

Individual cap
driver torque
control

Light touch
cap and decap
buttons

Long life
rechargeable
battery

Stand includes
convenient portrait
or landscape rack
placement



Low force dual
cap ejection
system

Cassettes compatible
with 24, 48 or 96 SBS
racks

Process screw
capped tubes in
row or columns